INSTRUCTIONS FOR USE VETSUTURE

p.2 NYLON
p.3 LENE
p.4 BOND
p.5 NEXT
p.6 PGA
p.7 fastPGA
p.8 PDX
p.9 PGC

NOTICES VETSUTURE

p.2 NYLON
p.3 LENE
p.4 BOND
p.5 NEXT
p.6 PGA
p.7 fastPGA
p.8 PDX
p.9 PGC

vetsuture® NYLON
NON-ABSORBABLE SURGICAL SUTURE (SYNTHETIC)
MONOFILAMENT (POLYAMIDE)

Description

Vetsuture® NYLON is a non-absorbable, sterile, surgical suture, prepared from long chain aliphatic polymers, Nylon 6 and Nylon 6.6. Vetsuture® NYLON is colored black with Logwood extract dye C 175290 conforming to U.S. Code of Federal Regulations. Vetsuture® NYLON is uncoated.

Indications

Vetsuture® NYLON is indicated for use in soft tissue approximation and or ligation, including use in Cardio-Vascular, Ophthalmic and Neurological procedures.

Actions

Vetsuture® NYLON, being a monofilament suture is remarkably smooth and easy to handle. It does not have interstices in its structure. It elicits an initial inflammatory reaction in tissues which is followed by gradual encapsulation of the suture by fibrous connective tissue. Vetsuture® NYLON is not absorbed, nor it is subjected to degradation or weakening by the action of tissue enzymes. Due to its relative biological inertness, it is recommended for use where the least possible suture reaction is desired. Because of its lack of adherence to tissue, Vetsuture® NYLON is effective as a pull out suture. Vetsuture® NYLON, being a monofilament has been successfully employed in surgical wounds which could subsequently become infected or contaminated where it can minimize later sinus formation and suture extrusion.

Warning

Do not re-sterilize. Discard open, unused sutures. Prolonged contact of this suture or any other suture with salt solutions such as those found in urinary and biliary tracts may result in calculus formation. Users should be familiar with surgical procedures and techniques involving non-absorbable sutures before employing Vetsuture® NYLON for wound closures, as risk of wound dehiscence may vary with the site of application and the suture material used. As any foreign material in the presence of bacterial contamination may enhance bacterial infectivity, acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.

Precautions

In handling this suture or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. As with any suture material, adequate knot security requires the accepted surgical techniques of flat and square ties, with additional throws as warranted by surgical circumstances and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting with Vetsuture® NYLON being a monofilament.

Adverse reactions

Adverse effects associated with the use of Vetsuture® NYLON include minimal initial inflammatory tissue reaction and transient local irritation at the wound site, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile juice etc. Vetsuture® NYLON may potentiate an existing infection.

Sterility

Vetsuture® NYLON suture are sterilized by ethylene oxide gas. Do not use if package is opened or damaged. Discard opened unused sutures.

Supply

Vetsuture® NYLON is available in several USP size in dozen boxes. Please visit WWW.VETSUTURE.COM to know our latest products and all NYLON sizes available.

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Do not re-use

Read instructions for use

Sterile unless the package is damaged or opened.
Method of sterilization : Ethylene oxyde

vetsuture® NYLON
SUTURE CHIRURGICALE NON-RESORBABLE (SYNTHETIQUE)
MONOFILAMENT (POLYAMIDE)

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**vetsuture® LENE**NON-ABSORBABLE SURGICAL SUTURE (SYNTHETIC)
MONOFILAMENT (POLYPROPYLENE)**Description**

vetsuture® NEXT is a non-absorbable synthetic monofilament suture made from polyvinylidene fluoride. vetsuture® LENE is a nonabsorbable, sterile, surgical suture, composed of Polypropylene, a synthetic linear Polyolefin. The molecular formula is (C3 H6)n. vetsuture® LENE is colored blue with Phthalocyanine Copper blue, Color Index No.74160, an U.S.F.D.A. approved dye. vetsuture® LENE is monofilament and it is uncoated. vetsuture® LENE meets all the requirements as per United States Pharmacopeia for nonabsorbable Surgical suture.

Indications

vetsuture® LENE is indicated for use in soft tissue approximation and or ligation, including use in plastic and ophthalmic procedures but not for use in cardiovascular and neurological tissues.

Actions

vetsuture® LENE being a monofilament suture is remarkably smooth and easy to handle. It does not have interstices in its structure. It elicits an initial inflammatory reaction in tissues which is followed by gradual encapsulation of the suture by fibrous connective tissue. vetsuture® LENE is neither absorbed, nor it is subjected to degradation or weakening by the action of tissue enzymes. Due to its relative biological inertness, it is recommended for use where the least possible suture reaction is desired. Because of its lack of adherence to tissue, vetsuture® LENE is effective as a pull out suture. vetsuture® LENE, being a monofilament has been successfully employed in surgical wounds which could subsequently become infected or contaminated where it can minimize later sinus formation and suture extrusion.

Warning

Do not resterilize. Discard open, unused sutures. Prolonged contact of this suture or any other suture with salt solution such as those found in urinary and biliary tracts, may result in calculus formation. Users should be familiar with surgical procedures and techniques involving nonabsorbable sutures before employing vetsuture® LENE for wound closures, as risk of wound dehiscence may vary with the site of application and the suture material used. As any foreign material in the presence of bacterial contamination may enhance bacterial infectivity, acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.

Precautions

In handling this suture or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holder. As with any suture material, adequate knot security requires the accepted surgical techniques of flat and square ties, with additional throws as warranted by surgical circumstances and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting with vetsuture® LENE being a monofilament.

Adverse reactions

Adverse effects associated with the use of vetsuture® LENE include minimal initial inflammatory tissue reaction and transient local irritation at the wound site, (calculus formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile juice etc. Like all foreign bodies, vetsuture® LENE may potentiate an existing infection.

Sterility

vetsuture® LENE suture are sterilized by ethylene oxide gas. Do not use if package is opened or damaged. Discard opened unused sutures.

Supply

vetsuture® LENE is available in several USP size in dozen boxes. Please visit WWW.VETSUTURE.COM to know our latest products and all LENE sizes available.

Storage

Recommandation storages conditions : below 25 Deg. centigrade away from moisture and direct heat. Do not use after expiry date.

Symbols used on the label

MFG Manufactured month & year

EXP Use until month & year

LOT Batch number

VET USE ONLY Veterinary use only

Do not re-use

Read instructions for use

STERILE unless the package is damaged or opened.
Method of sterilization : Ethilene oxyde

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**vetsuture® BOND**NON-ABSORBABLE SURGICAL SUTURE (SYNTHETIC)
BRAIDED AND COATED (POLYESTER)**Description**

vetsuture® BOND is a nonabsorbable, sterile, surgical suture, composed of fine filaments of Polyester. The fine polyester fibres are braided to produce firm suture that remains soft and pliable. vetsuture® BOND is colored green with D&C Green No.5, an U.S.F.D.A. approved dye. For added lubrication and smoothness, vetsuture® BOND is coated with highly biologically inert Polybutyrate or wax or with Teflon. vetsuture® BOND meets all the requirements, established by the United States Pharmacopoeia for Nonabsorbable Surgical suture.

Indications

vetsuture® BOND is indicated for use in soft tissue approximation and/or ligation, including use in plastic and ophthalmic procedures but not for use in cardiovascular and neurological tissues.

Actions

vetsuture® BOND being a braided suture, enables for secured knots. vetsuture® BOND is soft and pliable. It is very easy to handle. It elicits an initial inflammatory reaction in tissues which is followed by gradual encapsulation of the suture by fibrous connective tissue. vetsuture® BOND is not absorbed, nor it is subjected to degradation or weakening by the action of tissue enzymes. Due to its total biological inertness, it is recommended for use where the least possible suture reaction is desired. Because of its lack of adherence to tissue, vetsuture® BOND is effective as a pull out suture.

Warning

Do not resterilise. Discard open, unused sutures. Prolonged contact of this suture or any other suture with salt solutions such as those found in urinary and biliary tracts, may result in calculus formation. Users should be familiar with surgical procedures and techniques involving nonabsorbable sutures before employing vetsuture® BOND for wound closures, as risk of wound dehiscence may vary with the site of application and the suture material used. As any foreign material in the presence of bacterial contamination may enhance bacterial infectivity acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.

Precautions

In handling vetsuture® BOND suture or any other suture material care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. As with any suture material, adequate knot security requires the accepted surgical techniques of flat and square ties, with additional throws as warranted by surgical circumstances and the experience of the surgeon.

Adverse reactions

Adverse effects associated with the use of vetsuture® BOND include minimal initial inflammatory tissue reaction and transient local irritation at the wound site, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile juice etc. vetsuture® BOND may potentiate an existing infection.

Sterility

vetsuture® BOND sutures are sterilized by ethylene oxide gas. Do not use if package is opened or damaged. Discard opened unused sutures.

Supply

vetsuture® BOND is available in several USP size in dozen boxes. Please visit WWW.VETSUTURE.COM to know our latest products and all BOND sizes available.

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**vetsuture® NEXT**NON-ABSORBABLE SURGICAL SUTURE (SYNTHETIC)
MONOFILAMENT (POLYVINYLIDENE FLUORIDE)**Description**

vetsuture® NEXT is a non-absorbable synthetic monofilament suture made from polyvinylidene fluoride. vetsuture® NEXT meets all the general requirements established by the BP, the USP and the EP.

Indications

vetsuture® NEXT sutures are intended for use in procedures where non-absorbable monofilament sutures are indicated.

Actions

vetsuture® NEXT sutures exhibit excellent tissue tolerance, having virtually no tissue reaction.

Warnings

Users should be familiar with surgical procedures and techniques before employing vetsuture® NEXT. Acceptable surgical practice should be followed with respect to drainage and closure of contaminated or infected wounds.

Precautions

In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to the application of surgical instruments, such as forceps or needle holders. The smooth surface of the monofilament strand requires that knots be placed carefully and additional knots may be necessary to achieve knot security.

Sterility

vetsuture® NEXT suture are sterilized by ethylene oxide gas. Do not use if package is opened or damaged. Discard opened unused sutures.

Supply

vetsuture® NEXT is available in several USP size in dozen boxes. Please visit WWW.VETSUTURE.COM to know our latest products and all NEXT sizes available.

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Do not re-use

Read instructions for use

 STERILE Sterile unless the package is damaged or opened.
Method of sterilization : Ethilene oxyde**vetsuture®**
Sutures sold under trademark vetsuture®**vetsuture® NEXT**SUTURE CHIRURGICALE NON-RESORBABLE (SYNTHETIQUE)
MONOFILAMENT (FLUORIDE DE POLYVINYLIDENE)**Description**

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**vetsuture® PGA**ABSORBABLE SURGICAL SUTURE (SYNTHETIC)
BRAIDED AND COATED (POLYGLYCOLIC ACID)**Description**

PGA® is an absorbable, sterile, synthetic, surgical suture, composed of homopolymers of Glycolide (100%). PGA® is colored violet with D & C Violet No.2, conforming to U.S. Code of Federal regulations. For added lubrication and smoothness, PGA® is coated with unique combination of Polycaprolactone and Calcium Stearate. PGA® meets all the requirements, established by the United States Pharmacopeia for Absorbable Surgical suture.

Indications

PGA® is indicated for use in general soft tissue approximation and or ligation, including use in ophthalmic procedures but not for use in Cardio-Vascular and Neurological tissues.

Contra indications

PGA® suture being absorbable should not be used where extended approximation of tissue is required.

Actions

PGA® sutures get absorbed in the site of implantation, owing to the hydrolysis, where the polymer degrades to Glycolic acid. It is slowly absorbed and metabolized in the body. Absorption begins with loss of tensile strength followed by loss of mass. 80% of the original tensile strength is retained up to 14 days after surgery. 90% of the initial tensile strength is lost at the end of the fourth week post implantation. Absorption of PGA® suture is minimal upto 10% in two weeks, 25% of the mass gets absorbed in four weeks, upto 75% in two months and absorption is essentially complete between 75 and 90 days.

Warning

Do not re-sterilize. Discard open, unused sutures. Prolonged contact of this suture or any other suture with salt solutions such as those found in urinary and biliary tracts may result in calculus formation. As an absorbable suture, it may act transiently as a foreign body. As any foreign material in the presence of bacterial contamination may enhance bacterial infectivity, acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds. The use of this suture may be inappropriate in patients suffering from conditions which may delay wound healing.

Precautions

As this is an absorbable suture, the surgeon in the closure of sites undergoing expansion, stretching or distention, which may require additional support, should consider the use of supplemental non-absorbable sutures. Under some circumstances, notably orthopedic procedure, immobilization by external support may be employed at the discretion of the surgeon. Skin suture, which remain in place for more than 7 days may cause localized irritation and should be snipped off or removed. In handling PGA® suture or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. As with any suture material, adequate knot security requires the accepted surgical techniques of flat and square ties, with additional throws as warranted by surgical circumstances and the experience of the surgeon.

Adverse reactions

Adverse effects associated with the use of PGA® include allergic response in certain patients, transient local irritation at the wound site, transient inflammatory foreign body response, erythema and induration during the absorption process of subcuticular sutures.

Sterility

vetsuture® PGA sutures are sterilized by ethylene oxide gas. Do not use if package is opened or damaged. Discard opened unused sutures.

Supply

vetsuture® PGA is available in several USP sizes. The suture is supplied attached to various needle types in boxes of one dozen. PGA is also supplied in ligapack and cassette without needle.

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Symbols used on the label

MFG Manufactured month & year
 EXP Use until month & year
 LOT Batch number
 VET USE ONLY Veterinary use only

Do not re-use
 Read instructions for use
 STERILE Unless the package is damaged or opened.
 Method of sterilization : Ethylene oxide

Designed in FRANCE by
 oxybox vetsuture SARL
 54 rue Stendhal,
 75020, PARIS, FRANCE, EU
 www.vetsuture.com
 contact@vetsuture.com

**vetsuture®**

Sutures sold under trademark vetsuture®

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Copyright vetsuture 2010
 Date : 15 jan 2010
 Rev: 1 FR
 Doc : P-MO-PGA-1

vetsuture® fastPGA

ABSORBABLE SURGICAL SUTURE (SYNTHETIC)
BRAIDED AND COATED (POLYGLYCOLIC ACID)

Description

vetsuture® fastPGA is an absorbable, sterile, surgical suture, composed of homopolymers of Glycolide (100%). For added lubrication and smoothness, vetsuture® fastPGA is coated with unique combination of Polycaprolactone and Calcium Stearate. vetsuture® fastPGA meets all the requirements, established by the United States Pharmacopoeia for Absorbable Surgical suture. vetsuture® fastPGA is available as undyed.

Indications

vetsuture® fastPGA is intended for use in soft tissue approximation where only short-term wound support is required and where the rapid absorption would be beneficial. Due to the absorption profile vetsuture® fastPGA is useful for skin closure, particularly for difficult animals, surgery, episiotomies and closure of oral, mucosa, vetsuture® fastPGA is not intended for use in ophthalmic, surgical cardiovascular and neurological procedures.

Contra indications

vetsuture® fastPGA should not be used where extended approximation of tissues under stress is required or where wound support beyond 7 days is required due to the rapid loss of tensile strength.

Application

Sutures should be selected and implanted depending on patient condition, surgical experience, surgical technique and wound size. vetsuture® fastPGA typically falls off 7-10 days post operatively or can be wiped off successfully with sterile gauze. Normally the removal of the sutures is not required.

Actions

vetsuture® fastPGA sutures elicit a minimal to moderate initial inflammatory tissue reaction. Progressive loss of tensile strength and eventual absorption of vetsuture® fastPGA sutures occur by means of hydrolysis, where the polymer degrades to glycolic acid. It gets absorbed quickly and metabolized in the site of implantation. Absorption begins with loss of tensile strength followed by loss of mass. vetsuture® fastPGA shows that 5 days post implantation approximately 50% of the original tensile strength remains. All of the original tensile strength is lost by approximately 10-14 days post implantation. The absorption of vetsuture® fastPGA occurs thereafter and is essentially complete by 42 days.

Warning / Precautions / Interactions

Users should be familiar with surgical procedures and techniques involving absorbable sutures before employing vetsuture® fastPGA for wound closure as risk of wound dehiscence may vary with the site of application and the suture material used. Surgeons should consider the in vivo performance when selecting a suture. As with any foreign body, prolonged contact of this suture or any other suture with salt solutions such as those found in urinary and biliary tracts, may result in calculus formation. As an absorbable suture, it may act transiently as a foreign body. Acceptable surgical practice should be followed for the management of contaminated or infected wounds. As vetsuture® fastPGA is an absorbable suture material, the use of supplemental nonabsorbable sutures should be considered by the surgeon in the closure of sites which may undergo expansion, stretching or distention or which may require additional support. Skin sutures which must remain in place longer than 7 days may cause localized irritation and should be snipped off or removed as indicated. Consideration should be taken in the use of absorbable sutures in tissues with poor blood supply as suture extrusion and delayed absorption may occur. Subcuticular sutures should be placed as deeply as possible to minimize the erythema

Designed in FRANCE by
oxybox vetsuture S.A.R.L
54 rue Stendhal,
75020, PARIS, FRANCE, EU
www.vetsuture.com
contact@vetsuture.com

vetsuture®
Sutures sold under trademark vetsuture®

and induction normally associated with the absorption process. The use of vetsuture® fastPGA may be inappropriate in elderly, malnourished or debilitated patients, or patients suffering from conditions which may delay wound healing. In handling vetsuture® fastPGA suture or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. As with any suture material, adequate knot security requires the accepted surgical techniques of flat and square ties, with additional throws as warranted by surgical circumstances and the experience of the surgeon. Care should be taken to avoid damage when handling surgical needles. Grasp the needle in an area one third (1/3) to one half (1/2) of the distance from the attachment end to the point. Grasping in the point area could impair the penetration performance & cause fracture of the needle. Grasping at the butt or attachment end could cause bending or breakage of the needle. Reshaping needles may cause them to lose strength and make less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury. Discard used needles in Sharps containers.

Adverse reactions

Adverse reactions associated with the use of the vetsuture® fastPGA include: allergic response in certain patients, transient local irritation at the wound site, followed by moderate transient inflammatory foreign body response, erythema and induration during the absorption process of subcuticular sutures. Like all foreign bodies, vetsuture® fastPGA may enhance an existing infection.

Sterility

vetsuture® fastPGA sutures are sterilized by gamma radiation. Do not use if package is opened or damaged. Discard opened unused sutures.

Supply

vetsuture® fastPGA is available in several USP size in dozen boxes. Please visit WWW.VETSUTURE.COM to know our latest products and all fastPGA sizes available.

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Do not re-use

Read instructions for use

Sterile unless the package is damaged or opened.
Method of sterilization : gamma radiation

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Copyright vetsuture 2010
Date : 15 jan 2010
Rev 1 FR
Doc : P-MO-PFGA-1

**vetsuture® PDX**ABSORBABLE SURGICAL SUTURE (SYNTHETIC)
MONOFILAMENT (POLYDIOXANONE)**Description**

PDX® (Polydioxanone), monofilament synthétique absorbable suture is prepared from the polyester, poly (p-dioxanone). The empirical molecular formula of the polymer is (C₁₄H₂₀O₃)_n. Polydioxanone polymer has been found to be non antigenic, nonpyrogenic and elicits only a slight tissue reaction during absorption. It is colored violet, dyed with D&C Violet.

Indications

PDX® monofilament synthetic absorbable sutures are indicated for use in all types of soft tissue approximation, including use in cardiovascular tissue where growth is expected to occur and ophthalmic surgery. PDX® suture is not indicated in adult cardiovascular tissue, microsurgery and neural tissue. These sutures are particularly useful where the combination of an absorbable suture and extended wound support (up to six weeks) is desirable.

Contra indications

These sutures being absorbable, are not to be used where prolonged (beyond six weeks) approximation of tissues under stress is required and are not to be used in conjunction with prosthetic devices, i.e., heart valves or synthetic grafts.

Actions

Two important characteristics describe the in vivo performance of absorbable sutures: first, tensile strength retention and second, the absorption rate (loss of mass). PDX® synthetic absorbable suture has been formulated to minimize the variability of these characteristics and to provide wound support through an extended healing period. The results of implantation studies of PDX® monofilament suture in animals indicate that approximately 70% of its original strength remains two weeks after implantation. At four weeks post-implantation, approximately 50% of its original strength is retained and at six weeks, approximately 25% of the original strength is retained. Data obtained from implantation studies in rats show that the absorption of these sutures is minimal and occurs within 90th post-implantation day. Absorption is essentially complete within six months.

Warning

The safety and effectiveness of PDX® sutures have not been established in neural tissue, adult cardiovascular tissue or for use in microsurgery. Users should be familiar with the surgical procedures and techniques involving absorbable sutures before employing vetsuture® PDX® suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. Physicians should consider the in vivo performance when selecting a suture for use in patients. The use of this suture may be inappropriate in elderly, malnourished or debilitated patients or in patients suffering from conditions which may delay wound healing. Do not resterilize. Discard opened packages and unused sutures. As with any foreign body, prolonged contact of this suture or any other suture with salt solutions such as those found in urinary and biliary tracts, may result in calculus formation. As an absorbable suture, it may act transiently as a foreign body. Acceptable surgical practice should be followed for the management of contaminated or infected wounds. As this is an absorbable suture material the use of supplemental nonabsorbable sutures should be considered by the surgeon in the closure of the sites which may undergo expansion, stretching or distention or which may require additional support.

Designed in FRANCE by
oxybox vetsuture SARL
54 rue Stendhal,
75020, PARIS, FRANCE,
www.vetsuture.com
contact@vetsuture.com

**vetsuture®**
Sutures sold under trademark vetsuture®**Precautions**

Skin sutures which must remain in place longer than 7 days may cause localized irritation and should be snipped off or removed as indicated. Conjunctival and vaginal mucosal sutures remaining in place for extended periods may be associated with localized irritation and should be removed as indicated. Subcuticular sutures should be placed as deeply as possible to minimize the erythema and induration normally associated with absorption. Under circumstances notably orthopedic procedures, immobilization of joints by external support may be employed at the discretion of the surgeon. Consideration should be taken in the use of absorbable sutures in tissue with poor blood supply as suture extrusion and delayed absorption may occur. In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. vetsuture® PDX suture knots must be properly placed to be secure. Adequate knot security requires the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstances and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting monofilaments. Avoid prolonged exposure to elevated temperature. To avoid damaging needle points and suture areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling 'surgical needles' to avoid inadvertent needle sticks. Discard used needles in 'sharp' containers.

Adverse reactions

Due to prolonged suture absorption, some irritation and bleeding has been observed in the conjunctiva and mild irritation has been observed in the vaginal mucosa.

Sterility

vetsuture® PDX sutures are sterilized by ethylene oxide gas. Do not use if package is opened or damaged. Discard opened unused sutures.

Supply

vetsuture® PDX is available in several USP size in dozen boxes. Please visit WWW.VETSUTURE.COM to know our latest products and all PDX sizes available.

Storage

Recommended storage conditions : below 25 Deg. centigrade away from moisture and direct heat. Do not use after expiry date.

Symbols used on the label

MFG Manufactured month & year

EXP Use until month & year

LOT Batch number

Veterinary use only

Do not re-use

Read instructions for use

Sterile unless the package is damaged or opened.

Method of sterilization : Ethylene oxide

Copyright vetsuture 2010
Date : 15 Jan 2010
Rev: 1 EN
Doc : P-MO-PDX-1**vetsuture® PDX**SUTURE CHIRURGICALE RESORBABLE (SYNTHETIQUE)
MONOFILAMENT (POLYDIOXANONE)**Description**

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Designed in FRANCE by
oxybox vetsuture SARL
54 rue Stendhal,
75020, PARIS, FRANCE, EU
www.vetsuture.com
contact@vetsuture.com

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vetsuture® PGC (Poliglecaprone 25) is a monofilament synthetic absorbable suture, composed of a Polymer (Glycolide-co-caprolactone). vetsuture® PGC sutures are available undyed. Poliglecaprone 25 copolymer has been found to be nonantigenic, nonpyrogenic and elicits only a slight tissue reaction during absorption. vetsuture® PGC sutures comply with USP requirements, except for diameter.

Indications

vetsuture® PGC sutures are indicated for use in general soft tissue approximation, but not for use in cardiovascular or neurological tissues, microsurgery or ophthalmic surgery.

Contra indications

vetsuture® PGC sutures being absorbable, should not be used where extended approximation of tissue under stress is required, such as in fascia.

Actions

vetsuture® PGC suture is a monofilament which elicits a minimal acute inflammatory reaction in tissues and ingrowth of fibrous connective tissue. Progressive loss of tensile strength and eventual absorption of vetsuture® PGC sutures occurs by means of hydrolysis. Absorption begins as loss of tensile strength followed by loss of mass. Implantation studies in rats indicate that vetsuture® PGC suture retains approximately 50-60% of its original tensile strength 7 days post implantation and approximately 20-30% of its original tensile strength at 14 days post implantation. All of the original tensile strength is lost by 21 days post implantation. The absolute strength remaining 14 days post implantation meets or exceeds that historically observed with plain and chromic gut sutures. Absorption of vetsuture® PGC absorbable synthetic suture is essentially complete between 90 and 120 days.

Warning

Users should be familiar with the surgical procedures and techniques involving absorbable sutures before employing vetsuture® PGC suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. Physicians should consider the in vivo performance when selecting a suture for use in patients. The use of this suture may be inappropriate in elderly, malnourished or debilitated patients or in patients suffering from conditions which may delay wound healing. Do not resterilize. Discard opened packages and unused sutures. As with any foreign body, prolonged contact of this suture or any other suture with salt solutions such as those found in urinary and biliary tracts, may result in calculus formation. As an absorbable suture, it may act transiently as a foreign body. Acceptable surgical practice should be followed for the management of contaminated or infected wounds. As this is an absorbable suture material the use of supplemental nonabsorbable sutures should be considered by the surgeon in the closure of the sites which may undergo expansion, stretching or distention or which may require additional support.

Precautions

Skin sutures which must remain in place longer than 7 days may cause localized irritation and should be snipped off or removed as indicated. Subcuticular sutures should be placed as designed to minimize the erythema and induration normally associated with absorption. Under some circumstances, notably orthopedic procedures, immobilization of joints by external

suture may be employed at the discretion of the surgeon. Consideration should be taken in the use of absorbable sutures in tissue with poor blood supply as suture extrusion and delayed absorption may occur. In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. vetsuture® PGC suture knots must be properly placed to be secure. Adequate knot security requires the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstances and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting monofilaments. Avoid prolonged exposure to elevated temperature. To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in 'sharp' containers.

Adverse reactions

Adverse effects associated with the use of synthetic absorbable sutures include wound dehiscence, failure to provide adequate wound support in closure of the sites where expansion, stretching, or distention occur, failure to provide adequate wound support in elderly, malnourished or debilitated patients or in patients suffering from conditions which may delay wound healing, infection, localized irritation when skin sutures are left in place for more than 7 days, suture extrusion and delayed absorption in tissue with poor blood supply, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, and transitory local irritation at the wound site. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of blood-borne pathogens.

Sterility

vetsuture® PGC sutures are sterilized by ethylene oxide gas. Do not use if package is opened or damaged. Discard opened unused sutures.

Supply

vetsuture® PGC is available in several USP size in dozen boxes. Please visit WWW.VETSUTURE.COM to know our latest products and all PGC sizes available.

Storage

Recommended storage conditions : below 25 Deg. centigrade away from moisture and direct heat. Do not use after expiry date.

Symbols used on the label

MFG Manufactured month & year

EXP Use until month & year

LOT Batch number

VET USE ONLY Veterinary use only

Do not re-use

Read instructions for use

Sterile unless the package is damaged or opened.
Method of sterilization : Ethylene oxide

**vetsuture® PGC**SUTURE CHIRURGICALE RESORBABLE (SYNTHETIQUE)
MONOFILAMENT (POLIGLECAPRONE 25)**Description**

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